

Certificate of Conformity

Pharmatech SS Mounting Ring

Pharmatech confirms that the Mounting Ring enclosed has been manufactured from 316 Grade Stainless Steel, Satin Finish to Ra <1. Height 50mm. Solid GMP Rim. Fully Welded and Polished. It is Certified free from any imperfections:

Nominal Diameter: 250mm
Batch No: 08-250-02
Part No: 096/250
Material: 316 Stainless Steel

250 Ø

Approved by:

S Jennings

Date:

07.05.09

Pharmatech (a division of Terriva Limited)

2 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

www.pharmatech.co.uk email: info@pharmatech.co.uk