


# 1. CERTIFICATE OF CONFORMITY

2. Standards: **Code of Federal Regulations Title 21**
3. Equipment Description: **ø150mm ND Dust Cap Seal in Blue Silicone Rubber**
- 3.1 Pharmatech Part Reference: **095/150/B**
- 3.2 Batch Number: **13379/13125**
4. Manufacturer: **Connection Systems, 6 Forge Mills Park, Station Road Industrial Estate, Coleshill, Warwickshire. B46 1JH. U.K.**
5. This equipment and any acceptable variation thereto are specified in the **Technical Specification Datasheet**. This equipment is intended for use as a dust free connection of mounting rings and spigots with corresponding nominal diameter.
6. **Connection Systems** certifies that the goods supplied against the above purchase order have been designed and manufactured from **Silicone Rubber (PR110/50)** in line with the Technical Specification in a production process managed by a certified ISO 9001 Quality Assurance System. It has been inspected and is certified free from any imperfections.
7. Standards of manufacture are, but not limited to, the following:
- 7.1 **FDA CFR Title 21 Schedule 177.2600 - Rubber articles intended for repeated use.**
- 7.2 Tested in accordance with **USP Class VI** Biological testing.
- 7.3 Free from animal derived products or raw materials in line with EU Guidance **EMA/410/01**.
8. *It is the responsibility of the customer to ensure that this equipment is thoroughly cleaned before first use by their own internal validated cleaning process for Silicone products.*
9. This certificate relates only to the design and construction of the Equipment listed above and applies only for the intended use as specified in paragraph 5.

This is to certify that the equipment described herein ("Certified Equipment") has been manufactured to the specification identified above, in accordance with the standards and definitions dictated within Connection Systems Quality Assurance Manual. This certificate and test results obtained apply only to the certified equipment as submitted by Connection Systems.  
The certificate may not be used, in whole or in part, in any other document without Connection Systems prior written approval.

  
Shaun O'Riordan  
QA & Document Manager  
Date: 21 Oct 2019.

