

# Certificate of Conformity

## Pharmatech Dust Cap Seal

Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 250mm\*  
Batch No: 11/2996/4849  
Part No: 095/250/T  
Material: Silicone Rubber  
Colour: Translucent  
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:  
Title 21, Chapter 1, Schedule 177.2600

\*The Dust cap Seal is designed to be stretched over a mounting ring which is welded to a 250mm nominal bore pipe. Therefore the actual diameter of the seal will not be 250mm.

# 250 T

Approved by: .....



Date: 12th December 2011

**Pharmatech** (a division of Terriva Limited)

2 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)