

# Certificate of Conformity

## Connection Systems Dust Cap Seal

**Connection Systems** confirms that the Dust Cap Seal enclosed has been manufactured from CFR Compliant virgin, silicone rubber. This Dust Cap Seal has been inspected according to Connection Systems Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter:	100mm*
Batch No:	11885/11456
Part No:	995/100/B
Material:	Silicone Rubber
Colour:	Blue
Curing:	Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to CFR Standard:  
Title 21, Chapter 1, Schedule 177.2600  
U.S.P. Class VI

### TSE/BSE Statement:

We confirm that no materials derived from Animal Origin are used in the manufacture of these products , and no intermediates and/or auxiliary agents of animal origin have been used.

# 100 B

Approved by: ..... *Robinda* ..... Date: ..... *25 JAN 2018* .....

Name: ..... *SWAN O'RODAN* ..... Position: ..... *QA & DOCUMENTATION MANAGER* .....

### Connection Systems (a division of Terriva Limited)

Unit 6 Forge Mills Park  
Station Road Industrial Estate  
Coleshill, Warwickshire. B46 1JH  
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582  
[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)