

Certificate of Conformity

FDA IBC LID 'O' Ring Seal

Pharmatech confirms that the 'O' Ring Seal enclosed has been manufactured from FDA acceptable virgin silicone rubber. This 'O' Ring Seal has been inspected and is certified free from any imperfections.


Nominal Diameter: 585mm*
Tube Details: ID 5mm Wall Thickness 2.5mm
Batch No: 2961/4817
Pharmatech Part No: F004822 -1839
Material: Translucent Silicone Rubber

This 'O' ring Seal is certified as conforming to CFR Standard:
Title 21, Chapter 1, Schedule 177.2600

585 Ø

- The 'O' Ring Seal is designed for use with the lid for a Pharmatech IBC. The seal is made slightly undersize so it will stay in place on the lid when the lid is inverted.

IMPORTANT NOTICE: ensure that this 'O' Ring Seal is thoroughly cleaned before it is used for the first time.
Use your validated cleaning method for silicone products.

Approved by:  Date: 2nd December 2011
Name: Mr. M. J. Fenner Title: M.D.

Pharmatech (a division of Terriva Limited)
2 Forge Mills Park
Station Road Industrial Estate
Coleshill, Warwickshire. B46 1JH
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582
www.pharmatech.co.uk email: info@pharmatech.co.uk