

# Certificate of Conformity

## 'O' Ring Seal

### for 560mm diameter Drum Lid

Pharmatech confirms that the 'O' Ring Seal enclosed has been manufactured from FDA acceptable virgin silicone rubber. This 'O' Ring Seal has been inspected and is certified free from any imperfections.

Nominal Diameter: 560mm\*  
 Cord Diameter: 8.2mm  
 Batch No: 3096/5015  
 Pharmatech Part No: F000804  
 Colour: White

This 'O' ring Seal is certified as conforming to CFR Standard:  
 Title 21, Chapter 1, Schedule 177.2600

560 Ø

- The 'O' Ring Seal is designed for use with the lid for a Pharmatech 560mm ID drum. The seal is made slightly undersize so it will stay in place on the drum lid when the lid is inverted.

IMPORTANT NOTICE: ensure that this 'O' Ring Seal is thoroughly cleaned before it is used for the first time.  
 Use your validated cleaning method for silicone products.

Approved by: .....

*Col Beattie*

Date: .....

*3rd October 2013*

Name: .....

*Mr. Colin Beattie*

Title: .....

*Documentation Manager*

**Pharmatech** (a division of Terriva Limited)

4-5 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)