

Certificate of Conformity

Pharmatech Dust Cap Seal


Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 250mm*
Batch No: 3470/5698
Part No: 095/250/B
Material: Silicone Rubber
Colour: Blue
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:
Title 21, Chapter 1, Schedule 177.2600

*The Dust cap Seal is designed to be stretched over a mounting ring which is welded to a 250mm nominal bore pipe. Therefore the actual diameter of the seal will not be 250mm.

250 B

Approved by:  Date: 4th February 2013
Name: Mr. M. J. Fenner Position: MD

Pharmatech (a division of Terriva Limited)
4-5 Forge Mills Park
Station Road Industrial Estate
Coleshill, Warwickshire. B46 1JH
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582
www.pharmatech.co.uk email: info@pharmatech.co.uk