

# Certificate of Conformity

## Pharmatech Dust Cap Seal


Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter:	300mm*
Batch No:	4151/6245
Part No:	095/300/B
Material:	Silicone Rubber
Colour:	Blue
Curing:	Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:  
Title 21, Chapter 1, Schedule 177.2600

\*The Dust cap Seal is designed to be stretched over a mounting ring which is welded to a 300mm nominal bore pipe. Therefore the actual diameter of the seal will not be 300mm.

# 300 B

Approved by:  Date: 7th October 2013

Name: Mr. M. J. Fenner Position: MD

**Pharmatech** (a division of Terriva Limited)  
4-5 Forge Mills Park  
Station Road Industrial Estate  
Coleshill, Warwickshire. B46 1JH  
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582  
[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)