

Certificate of Conformity

Pharmatech Dust Cap Seal


Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 250mm*
Batch No: 5304/7505
Part No: 095/250/T
Material: Silicone Rubber
Colour: Translucent
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:
Title 21, Chapter 1, Schedule 177.2600

250 T

Approved by:



Date:

1st July 2015

Name:

DAVID EDWARDS

Position:

DOCUMENT CONTROL

Pharmatech (a division of Terriva Limited)
4-5 Forge Mills Park
Station Road Industrial Estate
Coleshill, Warwickshire. B46 1JH
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582
www.pharmatech.co.uk email: info@pharmatech.co.uk