

Certificate of Conformity

'O' Ring Seal

for 457mm diameter Drum Lid

Pharmatech confirms that the 'O' Ring Seals Fitted have been manufactured from FDA acceptable virgin silicone rubber. These 'O' Ring Seals have been inspected and are certified free from any imperfections.

Nominal Diameter: 457mm*
 Cord Diameter: 8mm
 Batch No: 5341/7457
 Pharmatech Part No: F000803

Colour: White

This 'O' ring Seal is certified as conforming to CFR Standard:
 Title 21, Chapter 1, Schedule 177.2600

We confirm that these items are not derived from animal origin, and during manufacture no intermediates and / or auxiliary agents which are of animal origin were used.

457 Ø

- The 'O' Ring Seal is designed for use with the lid for a Pharmatech 457mm ID drum. The seal is made slightly undersize so it will stay in place on the drum lid when the lid is inverted.

IMPORTANT NOTICE: ensure that this 'O' Ring Seal is thoroughly cleaned before it is used for the first time. Use your validated cleaning method for silicone products.

Approved by: Shaun O'Bordan Date: 17 Oct 16

Name: SHAUN O'BORDAN Position: Documentation Controller

Pharmatech (a division of Terriva Limited)

6 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

www.pharmatech.co.uk email: info@pharmatech.co.uk