

Certificate of Conformity

Pharmatech Dust Cap Seal

Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 200mm*
Batch No: 5752 / 9476
Part No: 095/200/T
Material: Silicone Rubber
Colour: Translucent
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:
Title 21, Chapter 1, Schedule 177.2600

TSE/BSE Statement:

We confirm that no materials derived from Animal Origin are used in the manufacture of these products , and no intermediates and/or auxiliary agents of animal origin have been used.

200 T

Approved by: *Michael* Date: 09 FEB 2018

Name: *Shawn O'Riordan* Position: QA & Document Manager

Pharmatech (a division of Terriva Limited)

6 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

www.pharmatech.co.uk email: info@pharmatech.co.uk