

Certificate of Conformity

Connection Systems Dust Cap Seal


Connection Systems confirms that the Dust Cap Seal enclosed has been manufactured from FDA Acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Connection Systems Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter:	250mm*
Batch No:	5995/8681
Part No:	995/250/T
Material:	Silicone Rubber
Colour:	Translucent
Curing:	Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:
Title 21, Chapter 1, Schedule 177.2600

250 T

Approved by:



Date:

13th January 2016

Name:

Mr. David Edwards

Position:

Document Control

Connection Systems (a division of Terriva Limited)

Unit 6 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

www.pharmatech.co.uk email: info@pharmatech.co.uk