

# Certificate of Conformity

## Connection Systems Dust Cap Seal


**Connection Systems** confirms that the Dust Cap Seal enclosed has been manufactured from FDA Acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Connection Systems Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 250mm\*  
Batch No: 6025/9338  
Part No: 995/250/B  
Material: Silicone Rubber  
Colour: Blue  
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:  
Title 21, Chapter 1, Schedule 177.2600

# 200 B

Approved by: .....



Date: .....

5th April 2016

Name: .....

Mr. David Edwards

Position: .....

Document Control

**Connection Systems** (a division of Terriva Limited)

Unit 6 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)