

# Certificate of Conformity

## 'O' Ring Seal

### for 750mm diameter Drum Lid

Pharmatech confirms that the 'O' Ring Seal enclosed has been manufactured from FDA acceptable virgin silicone rubber. This 'O' Ring Seal has been inspected and is certified free from any imperfections.

Nominal Diameter: 750mm\*  
Cord Diameter: 10mm  
Batch No: W0050879  
Pharmatech Part No: F000806  
Colour: White

This 'O' ring Seal is certified as conforming to CFR Standard:  
Title 21, Chapter 1, Schedule 177.2600

750 Ø

- The 'O' Ring Seal is designed for use with the lid for a Pharmatech 750mm ID drum. The seal is made slightly undersize so it will stay in place on the drum lid when the lid is inverted.

IMPORTANT NOTICE: ensure that this 'O' Ring Seal is thoroughly cleaned before it is used for the first time.  
Use your validated cleaning method for silicone products.

Approved by: .....

*S. Jennings*

Date: .....

*3rd. November 2009*

**Pharmatech** (a division of Terriva Limited)

2 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)