

# Certificate of Conformity

## 'O' Ring Seal

### for 457mm diameter Drum Lid

Pharmatech confirms that the 'O' Ring Seal enclosed has been manufactured from FDA acceptable virgin silicone rubber. This 'O' Ring Seal has been inspected and is certified free from any imperfections.

Nominal Diameter: 457mm\*  
Cord Diameter: 8mm  
Batch No: W004734  
Pharmatech Part No: F000803  
Colour: White

This 'O' ring Seal is certified as conforming to CFR Standard:  
Title 21, Chapter 1, Schedule 177.2600

457 Ø

- The 'O' Ring Seal is designed for use with the lid for a Pharmatech 400mm ID drum. The seal is made slightly undersize so it will stay in place on the drum lid when the lid is inverted.

**IMPORTANT NOTICE:** ensure that this 'O' Ring Seal is thoroughly cleaned before it is used for the first time.  
Use your validated cleaning method for silicone products.

Approved by: .....



Date: .....

02-09-09

**Pharmatech** (a division of Terriva Limited)

2 Forge Mills Park

Station Road Industrial Estate

Coleshill, Warwickshire. B46 1JH

Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582

[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@pharmatech.co.uk](mailto:info@pharmatech.co.uk)