

# **Certificate of Conformity**

## **Pharmatech Dust Cap Seal**

Pharmatech confirms that the Dust Cap Seal enclosed has been manufactured from FDA acceptable virgin, silicone rubber. This Dust Cap Seal has been inspected according to Pharmatech Standard Operating Procedure No: SOP 143 and is certified free from any imperfections:

Nominal Diameter: 400mm\*  
Batch No: W0050303  
Part No: 095/400/B  
Material: Silicone Rubber  
Colour: Blue  
Curing: Post cured for 4 hours at 200° C

This dust Cap Seal is certified as conforming to FDA Standard:  
Title 21, Chapter 1, Schedule 177.2600

\*The Dust cap Seal is designed to be stretched over a mounting ring which is welded to a 400mm nominal bore pipe. Therefore the actual diameter of the seal will not be 400mm.

# 400 B

Approved by: ..... *S. Jennings* ..... Date: ..... *29.04.09.* .....

**Pharmatech** (a division of Terriva Limited)  
2 Forge Mills Park  
Station Road Industrial Estate  
Coleshill, Warwickshire. B46 1JH  
Tel: +44 (0)1675 467 636 Fax: +44 (0)1675 467 582  
[www.pharmatech.co.uk](http://www.pharmatech.co.uk) email: [info@oharmatech.co.uk](mailto:info@oharmatech.co.uk)